

CI00402

- M1-M3
- RF Conducted **Immunity System**
- 100 Watts nominal
- 9kHz-400MHz

Features

Complete Testing Solutions to the following

- MIL-STD-461 CS114
- DO160 (Section 20) BCI Testing
- EN/IEC 61000-4-6
- IEC 60601-1-2
- EN 50130-4
- EN 61000-6-1/2
- EN 55024

The Model CI00402 a fully selfcontained state of the art system designed to test RF Conducted Immunity. The system contains a signal generator, spectrum



analyzer, 100W nominal AR amplifier 9 kHz to 400MHz, directional coupler, and emcware® control software. Everything is contained in a single housing, which eliminates setup issues. This system provides the versatility needed for

every test laboratory and equipment manufac-

Internal Test Specifications (See NOTE 1)
Cl00402
MIL-STD-461 (CS114)
DO160 (Sec 20 BCI Test)
IEC/EN 60601-1-2
IEC/EN 50130-4
IEC/EN 61326
IEC/EN 61000-6-1
IEC/EN 61000-6-2
IEC/EN 61000-4-6
CISPR 24/EN 55024
ISO 11452-4
Ford ES-XW7T-1A278-AC
GM GMW3097
BMW GS95002
Chrysler DC-11224
Renault 36-00-808
Peugeot B21 7110
NOTE 1: Specifications can be met using AR-

specified external accessories (injection probes, monitor probes, cal fixtures, CDN's, attenuators, etc.) Contact AR for further information.

turer. The RF amplifier, signal generator and spectrum analyzer can be used independently of the system. If special needs arise or standards were to change a larger amplifier can be connected to the system. The use of spectrum analyzer and signal generator may also be controlled by the software.

Option 1 adds data acquisition capability. Option 2 provides a laptop PC preloaded with the test software. Option 3 is a lower cost solution providing all the functions and features of the full Cl00402 base model without the cost of the internal amplifier.

The export classification for this equipment is EAR99. These commodities, technology or software are controlled for export in accordance with the U.S. Export Administration Regulations. Diversion contrary to U.S. law is prohibited.

The export classification for this equipment is EAR99. These commodities, technology or software are controlled for export in accordance with the U.S. Export Administration Regulations. Diversion contrary to U.S. law is prohibited.

Signal Generator Specifications		
Frequency range	9 kHz to 1.5 GHz	
resolution	0.01Hz	
Power range	-110 to +13 dBm	
resolution	0.01dB	
Modulation	AM, FM, Phase, Int Pulse,	
	Ext Pulse	

Spectrum Analyzer Specifications		
Frequency range	9 kHz to 1.5 GHz	
resolution	1 Hz	
RF Power CW (max)	20 dBm	
Atten = $30 dB$		
Resolution BW	10 Hz to 1 MHz	
Video BW	1 Hz to 3 MHz	
Amplitude Measure-	-110 dBm to +20 dBm	
ment Range	in 1 dB steps	
Preamplifier Gain	20 dB (nom)	
Sweep Time,	10 msec to 1500 sec	
span> 100 Hz		

AR RF/Microwave Instrumentation 160 School House Rd

Souderton, PA 18964 215-723-8181

For an applications engineer call:800.933.8181

www.arworld.us



Page 2

CI00402

- M1-M3
- RF Conducted Immunity System
- 100 Watts nominal
- 9kHz-400MHz

Specifications

RF Amplifier Specifications	
Frequency range	9 kHz to 400 MHz
Power rating	100 Watts Minimum
1dB compression	75 Watts Minimum
Harmonic Distortion	-20dBc at 75 Watts
Mismatch tolerance	100% of rated power without fold back. Will operate without damage or oscillation with any magnitude of source and load impedance.
Gain	50 dB minimum

Connections		
RF Out	Type N (front)	
Monitor Port In	Type N (front)	
Signal Generator Out	Type N (rear)	
RF Amp In	Type N (rear)	
Directional Coupler In	Type N (rear)	
RF Amp Out	Type N (rear)	
Pulse In	BNC (rear)	
Communication	USB B (rear)	
Directional Coupler Fwd Out	Type SMA (rear)	
Directional Coupler Fwd In	Type SMA (rear)	
Directional Coupler Rev Out	Type SMA (rear)	
Directional Coupler Rev In	Type SMA (rear)	

General		
Power	115/230 VAC, 50/60 Hz,	
	single phase 16A	
Breaker	2 pole, 20A	
Cooling	active cooling, air ventilation	
Environmental	10°C - 40°C	
conditions		
Dimensions,	50.3 x 53.3 x 55.1 cm	
	19.8 x 21.0 x 21.7 in	
Weight	49.9 kg (110 lb)	

Control Interface		
USB B	Connection to the computer	

Control Software		
emcware® Software Suite		
PC Requirements		
Computer	Intel Pentium 4, AMD Athlon 64 or better processor	
Operating system	Windows 7, Windows 8, Windows 10	
RAM	2 GB Minimum	
Screen Resolution	1024 x 768	
Ports	2 available USB 2.0 ports	
Software Require- ments	Microsoft Word/Excel 2007 or newer	

Options		
1	Data acquisition card	
2	Laptop PC with software preinstalled	
3	Amplifier removed; requires use of external amplifier.	

MODEL CONFIGURATIONS

MODEL	DESCRIPTION
CI00402M1	Includes Option 1
CI00402M2	Includes Option 2
CI00402M3	Includes Option 3

Accessory Kits

Application	Model	Description
IEC 61000-4-6 BCI Clamp 50Ω System	TK1000	Conducted immunity test kit containing all the attenuators, injection probes, calibration fixtures, calibration resistors, and termination resistors necessary for IEC testing.
IEC 61000-4-6 BCI Clamp 150Ω System	TK1001	Conducted immunity test kit containing all the attenuators, injection probes, calibration fixtures, calibration resistors, and termination resistors necessary for IEC testing.
IEC 61000-4-6 EM Clamp 50Ω System	TK1002	Conducted immunity test kit containing all the attenuators, injection probes, calibration fixtures, calibration resistors, and termination resistors necessary for IEC testing.
MIL-STD-461 DO160	TK2000	Conducted immunity test kit containing all the attenuators, injection probes, calibration fixtures, calibration resistors, and termination resistors necessary for MIL-STD-461 and DO160 testing.
Automotive	TK3000	Conducted immunity test kit containing all the attenuators, injection probes, calibration fixtures, calibration resistors, and termination resistors necessary for automotive testing.